AUSTRIA

Bauerfeind Ges.m.b.H.
Hainburger Straße 33
1030 Wien
P +43 (0) 800 4430-130
F +43 (0) 800 4430-131
E info@bauerfeind.at

BENELUX

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en english

These instructions for use can also be downloaded from www.bauerfeind.com/downloads.

Dear customer,

Thank you for choosing a product from Bauerfeind. The MalleoTrain® S open heel you have purchased is a high-quality medical product.

Please read through all of these **instructions for use** carefully, because they contain important information about the use, indications, contraindications, application, care and disposal of the MalleoTrain® S open heel. Please keep these instructions safe, as you may want to read them again later. If you have any questions, please contact your physician or medical retailer.

Intended use

The MalleoTrain® S open heel is an active support for the ankle. This support with an open heel will enable you to tread on the floor with a feeling of stability and better grip. It has been specifically developed for sports in which bare feet come into contact with the floor, such as martial arts, gymnastics or dancing. In addition to the support, the strap enables users to take advantage of the modes of action offered by sports tape. A strap helps to improve the way wearers feel the movement of their foot. Like sports

tape, this strap helps to optimize the way the body's center of gravity is guided via the foot. This helps to prevent injuries and speed up the healing process.

Indications

- Ligament insufficiency
- Post-operative rehabilitation
- Slight sprains
- Relapse prevention, particularly during sports activities
- Prevention

Risks of using this product

⚠ Caution*

- MalleoTrain® S open heel must only be worn in accordance with the specifications contained in these instructions and used for the areas of application (purpose, places of use) listed.
- To ensure the optimum fit of the MalleoTrain® S open heel, the size of MalleoTrain® S open heel required must be determined by a trained specialist¹.
- The sewn-in label of the MalleoTrain® S open heel, which contains information on the product name, size, manufacturer, washing instructions and CE mark, is located on the inside of the knitted fabric. Only a trained specialist¹ may carry out the (initial) fitting of the support and provide instruction.
- Use in combination with other products, e.g. as part of compression treatment (compression stockings), must be discussed beforehand with your physician. The benefits of this product are felt primarily during physical activity. Take off the MalleoTrain® S open heel during extended rest phases (e.g. sleeping, long periods of sitting down).
- Should you notice any unusual changes (e.g. increased symptoms), please contact your physician without delay.
- No product liability is accepted in the event of improper use. Do not allow the product to come into contact with ointments, lotions or substances containing grease or acids.
- Side effects that affect the entire body have not been reported to date. Correct fitting is assumed. Do not put the MalleoTrain® S open heel on too tightly as this could lead to local pressure symptoms. In rare cases, constriction of blood vessels and nerves may occur. If this happens, check the fit of your support, loosen the strap and have the size of your MalleoTrain® S open heel checked if necessary.

 Contains natural rubber (latex (Elastodiene)) which can cause skin irritation in persons with an allergy. Take off the support if you notice any allergic reaction.

Note: Never expose the MalleoTrain® S open heel to direct heat (e.g. heaters, sunlight, in your car)! This may cause damage to the material, which can impair the effectiveness of the MalleoTrain® S open heel.

Contraindications

No clinically significant hypersensitive reactions have been reported to date. If you have any of the following conditions such aids should only be fitted and applied after consultation with your physician:

- Skin disorders/injuries in the relevant part of the body, particularly if inflammation is present. Likewise, any raised scars with swelling, redness and excessive heat build-up
- Varicose veins
- Impaired sensation and circulatory disorders of the legs (e.g. diabetes mellitus)
- Impaired lymph drainage – including soft tissue swellings of uncertain origin located away from the fitted support

Application instructions

Putting on the support

- 1 Pull on the support so that the start of the strap is located on the inside of the foot. The start of the strap defines the direction in which the strap needs to be wound around the foot. When applying the strap, make sure that it sits without folds.

- 2 Bring the strap under the sole of your foot toward the outside of your foot. By pulling the strap you can stabilize the foot and ankle comfortably. To ensure the optimum fit of the strap, your lower leg must be perpendicular (at a right angle) to the foot during all the following steps. 3 Now bring the strap diagonally across the back of the foot and back around ankle.

- 4 5 Guide the strap diagonally forward again to the inside of the foot.

- 6 Pull the strap under the sole of your foot toward the outside of the foot and over the back. 7 Now secure the end of the strap by fastening it to the section of the strap underneath it. Alternatively, you can also use the indications on the support to guide you. Blue dots on the titanium-colored knitted fabric indicate where the strap is supposed to go. When wrapping the strap around your foot, follow the dots knitted into the fabric, starting with one dot, going to two dots and then to three dots.

Taking off the support

Undo the strap. To avoid damage to the support, fasten the end of the strap to another part of the strap. Take hold of the upper edge of the support and turn it down.

Cleaning instructions

Fasten the Velcro fastening to avoid damage. Wash the support separately in a lingerie bag with a mild detergent at 30 °C. Wash the MalleoTrain® S open heel regularly so that the knitted fabric retains its compressive properties.

Notes on reuse

The MalleoTrain® S open heel is intended for your treatment only and has been specifically selected according to your measurements. For this reason, do not pass the MalleoTrain® S open heel on to other people.

Warranty

The legal provisions of the country in which the product was purchased apply. Country-specific warranty regulations that govern the relationship between the retailer and the buyer also apply where relevant. If you believe that a claim may be made under the warranty, please first contact the retailer from which you purchased the product directly.

Please do not make any changes to the product yourself. This particularly applies to individual adjustments made by specialists. Please follow our instructions for use and care. These are based on our experience gained over many years and ensure the functionality of our medically effective products for a long time. Only aids that function optimally can help you. Failure to follow the instructions may also restrict the warranty.

Disposal

You can dispose of the product in accordance with national legal requirements.

Version: 2017-06

*Note on the risk of personal injury (risk of injury, harm to health and accidents) or damage to property (damage to the product).

¹ A specialist is any person who is authorized according to the state regulations for fitting and instruction in the use of supports which are relevant to you.

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