**AUSTRIA**

Bauerfeind Ges.m.b.H.  
Hainburger Straße 33  
1030 Wien  
P +43 (0) 800 4430-130  
F +43 (0) 800 4430-131  
E info@bauerfeind.at

**BENELUX**

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**These instructions for use can also be downloaded from:**  
[www.bauerfeind.com/downloads](http://www.bauerfeind.com/downloads)

## Dear Customer,

Thank you for choosing a product from Bauerfeind. The SecuTec OA knee orthosis<sup>1</sup> is a high-quality medical and therapeutic product. Please read the **instructions for use** carefully and if you have any questions contact your physician, medical retailer, or our Technical Service department.

### Intended purpose

The SecuTec OA orthosis is an approved therapeutic aid for supporting your treatment. It secures and provides relief for your knee during therapy and can be used in conjunction with physiotherapy to help reduce pain. By correcting the leg alignment, targeted relief for the knee joint and its ligaments can be achieved for the relevant indication.

**The SecuTec OA must be individually adjusted for the patient by a trained specialist<sup>2</sup> in order to ensure maximum effectiveness and optimum wearing comfort.**

### Indications

- Medial or lateral osteoarthritis of the knee
- Postoperative stabilization and relief

### Risks of using this product

#### Important information

- The SecuTec OA orthosis is a prescribed product that must be worn under a physician's guidance. The product must be worn in accordance with these instructions for use and for the listed indications.
- No product liability is accepted if the product is not used or adjusted correctly.
- No improper modifications may be made to the product. Failure to comply with this requirement may adversely affect product performance, thereby excluding any product liability.
- Should symptoms increase or should you notice any unusual changes, please contact your physician without delay.
- Combining the orthosis with other products (e.g. compression stockings) must be discussed beforehand with your physician.
- Do not allow the product to come into contact with ointments, lotions, or substances containing grease or acids.
- If the orthosis is used in salt water or chlorinated water, it should be rinsed with plain water afterwards.
- Should any problems with the fit of the orthosis arise, for instance if the shape of your leg changes (e.g. due to an increase or reduction in muscle), or should the product have a functional defect, please contact your medical supply retailer immediately.
- Side effects involving the body as a whole have not been reported to date. Correct usage / fitting is assumed. Any medical aids – supports and orthoses\* – applied externally to the body can, if tightened excessively, lead to local pressure symptoms or, rarely, constrict blood vessels or nerves.

### Contraindications

No clinically significant hypersensitivity reactions have been reported to date. If you have any of the following conditions, the product should only be fitted and worn after consultation with your physician:

1. Skin disorders / injuries in the relevant part of the body, particularly if inflammation is present. Likewise, any raised scars with swelling, redness, and excessive heat build-up.
2. Varicose veins (varicosis).
3. Impaired sensation and circulatory disorders of the legs / feet (e.g. diabetes mellitus)
4. Impaired lymph drainage – including soft tissue swellings of uncertain origin located away from the fitted aid.

### Application instructions

#### Size calculation

The orthosis is selected and fitted by a specialist after measuring the patient. You can refer to the size chart on the packaging for the correct size.

#### Procedure for the first fitting of the SecuTec OA knee orthosis

**Adjustment options for the hinge (to be performed by a trained specialist only)**

The hinge cover can be opened via a screw fastening. This cover fastens the movement-limiting wedges (extension – front / flexion – rear) and also protects the hinge mechanism against dirt.

Wedges are pre-installed at the factory. They primarily provide fixed stop absorption during extension and they do not limit the range of motion.

- Opening the screw fastening (quarter turn until the arrow points at the "open" symbol ) allows the hinge protection to be removed and the desired limiting wedges to be changed (1 – a – c).

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- The desired extension and flexion wedges are selected in pairs from the enclosed frame and are inserted into the orthosis hinges on both sides. The following limitations are possible:  
Extension: 0°, 10°, 20°, 30°, 45°; flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°; fixation: 0°, 10°, 20°, 30°, 45°.
- Make sure that the screw fastening is showing the  symbol and click the cover back on to the hinge. When the wedges are flush with the surface of the cover, the screw fastening can be turned back to .

### Adjustment (only by a trained specialist)

#### Preparation:

- The orthosis has special gel cushions on one side in order to ensure optimal pressure distribution.  
At the factory, the gel cushions for **correcting a varus knee (bandy legs)** or **medial osteoarthritis of the knee** are attached to the lateral hinge of the orthosis.  
In order to **correct a valgus knee (knock knees)** or **lateral osteoarthritis of the knee**, the gel cushions need to be exchanged.
  - The oblong gel cushion on the lateral side of the orthosis frame below the knee has a cut-out for the head of the fibula. Make sure that this is in the correct position.
  - Mark the compromise pivot axis of the knee joint (as defined by Nietert) on the leg.
  - The pivot point of the orthosis hinge lies in the center of the screw fastening on the hinge cover.
  - When fitting the orthosis for the first time, start by unclicking all four quick-fit fasteners and loosening the two other straps.
  - With the knee in a slightly flexed, unloaded position, put on the SecuTec OA in such a way that the pivot point of the orthosis hinge corresponds to the compromise pivot axis (as defined by Nietert).
  - Check the fit of the shell parts when they are laid flat to make sure their circumference and longitudinal direction are correct and reshape them if necessary. A bending iron with protective caps should be used for any extensive bending work.
  - The calf shell is fitted to the anatomy of the lower leg by adjusting the slant.
  - The position of the strap cushions can be altered as required.
  - Now close the quick-fit fasteners that are close to the hinge (#1 and #2), followed by the ones further away (#3 and #4) , and tighten the straps using the Velcro fastening .
  - Check once again that the orthosis hinges are in the correct position and tighten the Velcro straps close to the hinge (without a quick-fit fastener).
  - Adjust the hinges on the lower leg frame to the shape of the leg; a slight correction position may be set here if required .
  - A 3 mm Allen key can be used for a valgus / varus setting of up to 25° in each case. The direction of rotation is marked on the orthosis .
  - Please note: Extreme settings cause the patient pain and often offer no therapeutic benefit.
  - Turn the screw on the hinge in the direction in which the lower part of the frame should move. A scale and a groove in the screw will help to guide you here. The distance between two marks corresponds to a correction of 2°.
  - Make the same setting on the other side in order to ensure that the shift is parallel; this results in an optimal application of force and ensures that the hinges are parallel  (varus)  (valgus).
  - The gel condylar cushion can also be adjusted to the anatomy of the knee by using insertable wedges to increase the size of the contact surface for the correction. To make this adjustment, insert the wedges supplied (3 and / or 6 mm) from the front until the cushion is lying flat against the knee .
- ### Putting on the SecuTec OA
- Since the straps and their length have been adjusted by the specialist, these do not necessarily need to be undone. To put on the orthosis, just undo the four quick-fit fasteners.
  - Put on the orthosis by holding it in front of you, placing your foot down through the frame, and then pulling the frame up to your knee .
  - Bend your leg slightly (approx. 30°). With your fingers, feel for the position of the kneecap. The center of the orthosis hinges is positioned at approximately the same height as the center of the kneecap.
  - The grips of the quick-fit fasteners also have numbers (1 – 4) imprinted underneath and these indicate the best order for doing up the straps .
  - Click all of the quick-fit fasteners on the straps into the eyelets, starting at the front with the fastener on the lower leg close to the hinge (#1), followed by the one on the thigh close to the hinge (#2). Next, close the fasteners further away from the hinge, starting with the one on the lower leg first (#3), followed by the one on the thigh (#4).
  - If the straps are too loose or too tight, you can tighten or loosen them slightly at this point.  
Attention: Straps that are too loose may cause the orthosis to slip.
  - After all the straps have been closed, the orthosis must lie evenly when the muscles are tensed.
  - Once the orthosis has been put on, check once more that it is in the correct position .

### Taking off the SecuTec OA

- Start by taking the quick-fit fasteners out of the eyelets on the back of the thigh frame of the orthosis.
- Next, undo the quick-fit fasteners of the straps at the front on the lower leg part.
- The front top strap and the rear bottom strap can remain closed.
- Once you have undone all of the quick-fit fasteners, push the orthosis down toward your foot and step out of the frame.

### Cleaning instructions

**Never expose the orthosis to direct heat (e.g. heater, sunlight, etc.)!**

The aluminum parts of the orthosis are coated and can be cleaned using pH-neutral soap.

We recommend washing the cushions and straps by hand at 30 °C.

This product has been tested by our own quality control system.

However, should you have any complaints, please contact your local medical supply retailer.

### Where to apply the orthosis

In accordance with the indications (knee). Please see intended purpose.

### Maintenance instructions

If handled and cared for correctly, the product is practically maintenance-free.

### Assembly and fitting instructions

The SecuTec OA is supplied as a pre-assembled orthosis in standard sizes. The first fitting, which includes the hinge adjustments and the static fitting, must be carried out by a trained specialist.

### Technical specifications / parameters

The SecuTec OA is a functional orthosis for the knee joint. It consists of a frame, hinges, straps, quick-fit fasteners, and cushions.

### Notes on reuse

The SecuTec OA is intended for the treatment of one patient only.

### Disposal

- This product can be disposed of according to national laws.
- For information regarding the options available for returning the packaging to ensure it is correctly recycled, please observe the instructions from the sales outlet where you purchased this product.

Version: 2018-06

<sup>1</sup> Orthosis = orthopedic aid used to stabilize, relieve, immobilize, control, or correct the limbs or the torso

<sup>2</sup> A specialist is any person who is authorized according to the state regulations for fitting and instruction in the use of supports that are relevant to you.

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