SONDES DOPPLER

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
Product

Waterproof vascular probe is used for monitoring arteries and veins.

Optional: 5.0 MHz/8.0 MHz

Ultrasound

Nominal frequency: 5.0MHz/8.0 MHz

Working frequency: 5.0MHz/8.0 MHz $\pm 10\%$

$P_{\text{max}} < 1\text{MPa}$

$I_{\text{ob}} < 20\text{ mW/cm}^2$

$I_{\text{sata}} < 10\text{ mW/cm}^2$

$I_{\text{spta}} < 100\text{ mW/cm}^2$

Effective area for transmitter:

$32\text{ mm}^2 \pm 15\%$ (4 MHz / 5.0 MHz)

$14\text{ mm}^2 \pm 15\%$ (8.0 MHz)

Degree of protection against harmful ingress of water

Protected against the effects of continuous immersion in water with 1m depth for 24 hours.

Degree of safety in presence of flammable gases

Equipment not suitable for use in presence of flammable gases.

Working system: Continuous running equipment.

Working mode: Continuous wave Doppler.

General operation

1. Connect an optional probe to a vascular Doppler main unit.

2. Apply a liberal amount of gel on the site to be examined (see Figure 1 Probe sites). Place the probe at 45° to the skin surface over the vessel to be examined.

3. Adjust the position of the probe to obtain the loudest audio signal. High pitched pulsatile sounds are emitted from arteries while veins emit a non-pulsatile sound similar to a rushing wind.

WARNING: please use the vascular Doppler main unit developed by our company or the probes may be damaged.

Figure 1 Probe sites

Cleaning and disinfecting

1. Wipe the probe with a soft cloth to remove any remaining coupling gel. Clean with soap and water only.

2. After cleaning may use 70 % ethanol to disinfect. Wipe the probe with a clean, dry cloth to remove any remaining moisture.

WARNING

1. Do not use this device in or around the eye.

2. Don't use strong solvent, for example, acetone.

3. Never use an abrasive such as steel wool or metal polish.

4. This device is not explosion-proof and can not be used in the presence of flammable anaesthetics.

5. After the service life, please return the products to the manufacture or dispose the products according to local regulations.

6. Do not use high temperature sterilizing process and E-beam or gamma radiation sterilization.

7. This device is not intended for treatment. The intended use is for clinical use. If the result is distrustful, please use other methods to verify immediately.

Figure 1 Probe Sites

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
*SPENGLER,
THE QUALITY BRAND FOR HEALTH PROFESSIONALS*