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SCOOTER S19

HEARTWAY MEDICAL PRODUCTS Co., Ltd.

TAIWAN

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SOMMAIRE

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CHARGING INSTRUCTIONS

To recharge the batteries, follow the steps below:

- ⇒ Place your power scooter close to a standard electrical wall outlet.
- ⇒ Remove your key to turn the power OFF
- ⇒ Slide the charger port door open.

- ⇒ Plug the XLR connector of the charger to the charger port.
- ⇒ Plug the other end of power cord into a standard wall outlet.
- ⇒ When charging is completed, battery capacity indicator is shown.
- ⇒ Disconnect the charger power cord from the wall outlet when the batteries are fully charged.

Recharge battery only when the key is in off position. When indicator is in low status, this confirms the battery needs recharging.

Note:

- ⇒ Always charge your batteries in well ventilated areas.
- ⇒ The charger is intended for indoor use only. Protect from moisture.
- ⇒ For maximum performance, it is recommended that you replace both batteries at the same time if the batteries are weak.
- ⇒ If the vehicle will not be used for a long period of time, arrange to have the batteries recharged at least once every month to avoid deterioration of the batteries.

According to the battery type and condition of the batteries, they usually can be fully charged in 4-10 hours. This will be indicated when the status light in the battery charger side panel turns green. Charging the battery longer than necessary will not harm the battery. We recommend that you charge the batteries for 8 to 10 hours after daily use.

BATTERY INDICATOR

A battery indicator is to provide information about the travelling range remaining.

Indicator	Travelling Range Remaining
Full	15 km Traveling Range Remaining
Green Area	7~15 km Travelling Range Remaining
Yellow Area	4~7 km Travelling Range Remaining
Red Area	1~4 km Travelling Range Remaining
End	0 Kilometer

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HEARTWAY MEDICAL PRODUCTS CO.,

Ltd.

Explanatory sheet about safety of product for transportation
(Safety Data Sheet for transportation)

1. Basic item

Product name	Lithium ion cell
Product code	Refer to Table 1.
Product model name	Refer to Table 1.
Manufacturer	SANYO Electric Co., Ltd. Energy Company
Address	222-1 Kaminaizen Sumoto Hyogo, Japan
Department in charge	Lithium-Ion Battery Business Unit Battery System Management Department Technical Administration Group
Phone number	+81-799-23-3931
E-mail	joho_gijutsu@gg.jp.panasonic.com


2. Product information

Basic composition of the product

This product is a Lithium ion (polymer) single with or without tube casing.

3. Safety Information

1. The Watt-hour rating of these Lithium ion (polymer) cells is not more than 20Wh.
Refer to Appendix "SDS (SDS-IAT-00002)".
2. These Lithium ion (polymer) cells can be treated as non-dangerous goods, this means they are not categorized to dangerous goods defined by UN3480, under the UN Recommendation on the Transport of Dangerous Good, Special Provision 188.
3. SANYO certifies the Lithium ion cells have passed and satisfied the UN Manual of Tests and Criteria Part III, sub-section 38.3 testing in Sanyo shipping.
4. Packages of these Lithium ion (polymer) cells satisfy the following conditions when SANYO ships.
 - (1) The product name "Lithium ion batteries" and how to deal with the damage of the package are written on the label and attached document.
 - (2) The packages have passed the drop test from the height of 1.2m.
 - (3) Package weight is not more than 10kg.
 - (4) IATA DGR 53rd edition Packing Instruction 965.
5. Since the package weight of these lithium ion (polymer) cells is not more than 10kg in Sanyo shipping, these lithium ion (polymer) cells are exempt from IATA DGR 53rd edition Special Provision A99.
6. Since these Lithium ion (polymer) cells have passed and satisfied the UN Manual of Tests and Criteria Part III, sub-section 38.3 testing, these Lithium ion (polymer) cells are exempt from IATA DGR 53rd edition Special Provision A88.

General Manager H. Kuroda
Technical Administration Group
Battery System Management Department
Lithium-Ion Battery Business Unit
SANYO Electric Co., Ltd. Energy Company

Safety data sheet for product

1. PRODUCT AND COMPANY IDENTIFICATION

- Product name: Lithium ion battery cell
- Product code: None (All models Sanyo manufactured and whose capacity is less than or equal to 4.95Ah)
- Company name: Sanyo Electric Co., Ltd. Energy Company
- Address: 222-1, Kaminaizen, Sumoto City, Hyogo, Japan
- Telephone number: +81-799-24-4111
- Fax number: +81-799-23-2879
- Emergency telephone number: [Weekday] +81-799-23-3931
[Night and holiday] +81-799-24-4131

2. COMPOSITION / INFORMATION ON INGREDIENTS

- Substance or preparation: Preparation
- Information about the chemical nature of product: *1

Common chemical name / General name	CAS number	Classification and hazard labeling
Lithium transition metal oxide ($\text{Li}[\text{M}]_m[\text{O}]_n \cdot 2$)	12190-79-3 12057-17-9 182442-95-1	
Iron	7439-89-6	
Aluminum	7429-90-5	
Graphite (Natural graphite) (Artificial graphite)	7782-42-5 7440-44-0	
Copper	7440-50-8	
Organic electrolyte	-	Inflammable liquid

*1 Not every product includes all of these materials.

*2 The letter M means transition metal and candidates of M are Co, Mn, Ni and Al. One compound includes one or more of these metals and one product includes one or more of the compounds. The letter m and n means the number of atoms.

3. HAZARDS IDENTIFICATION

For the battery cell, chemical materials are stored in a hermetically sealed metal case, designed to withstand temperatures and pressures encountered during normal use. As a result, during normal use, there is no physical danger of ignition or explosion and chemical danger of hazardous materials' leakage.

However, if exposed to a fire, added mechanical shocks, decomposed, added electric stress by miss-use, the gas release vent will be operated. The battery cell case will be breached at the extreme, hazardous materials may be released.

Moreover, if heated strongly by the surrounding fire, acrid gas may be emitted.

- Most important hazard and effects

Human health effects:

Inhalation: The steam of the electrolyte has an anesthesia action and stimulates a respiratory tract.

Skin contact: The steam of the electrolyte stimulates a skin. The electrolyte skin contact causes a sore and stimulation on the skin.

Eye contact: The steam of the electrolyte stimulates eyes. The electrolyte eye contact causes a sore and stimulation on the eye. Especially, substance that causes a strong inflammation of the eyes is contained.

Environmental effects: Since a battery cell remains in the environment, do not throw out it into the environment.

- Specific hazards:

If the electrolyte contacts with water, it will generate detrimental hydrogen fluoride.

Since the leaked electrolyte is inflammable liquid, do not bring close to fire.

4. FIRST-AID MEASURES

Spilled internal cell materials

- Inhalation:
Make the victim blow his/her nose, gargle. Seek medical attention if necessary.
- Skin contact:
Remove contaminated clothes and shoes immediately. Wash extraneous matter or contact region with soap and plenty of water immediately.
- Eye contact:
Do not rub one's eyes. Immediately flush eyes with water continuously for at least 15 minutes. Seek medical attention immediately.

A battery cell and spilled internal cell materials

- Ingestion:
Make the victim vomit. When it is impossible or the feeling is not well after vomiting, seek medical attention.

5. FIRE-FIGHTING MEASURE

- Suitable extinguishing media: Plenty of water, carbon dioxide gas, nitrogen gas, chemical powder fire extinguishing medium and fire foam.
- Specific hazards: Corrosive gas may be emitted during fire.
- Specific methods of fire-fighting: When the battery burns with other combustibles simultaneously, take fire-extinguishing method which correspond to the combustibles. Extinguish a fire from the windward as much as possible.
- Special protective equipment for firefighters:
 - Respiratory protection: Respiratory equipment of a gas cylinder style or protection-against-dust mask
 - Hand protection: Protective gloves
 - Eye protection: Goggle or protective glasses designed to protect against liquid splashes
 - Skin and body protection: Protective cloth

6. ACCIDENTAL RELEASE MEASURES

Spilled internal cell materials, such as electrolyte leaked from a battery cell, are carefully dealt with according to the followings.

- Precautions for human body:
Remove spilled materials with protective equipment (protective glasses and protective gloves). Do not inhale the gas as much as possible. Moreover, avoid touching with as much as possible.
- Environmental precautions: Do not throw out into the environment.
- Method of cleaning up: The spilled solids are put into a container. The leaked place is wiped off with dry cloth.
- Prevention of secondary hazards: Avoid re-scattering. Do not bring the collected materials close to fire.

7. HANDLING AND STORAGE

- Handling
 - Technical measures:
 - Prevention of user exposure: Not necessary under normal use.
 - Prevention of fire and explosion: Not necessary under normal use.
 - Precaution for safe handling: Do not damage or remove the external tube.
 - Specific safe handling advice: Never throw out cells in a fire or expose to high temperatures. Do not soak cells in water or seawater. Do not expose to strong oxidizers. Do not give a strong mechanical shock or fling. Never disassemble, modify or deform. Do not connect the positive terminal to the negative terminal with electrically conductive material. Do not use in the place temperature can be high. In the case of charging, use only dedicated charger or charge according to the conditions specified by Sanyo.
- Storage
 - Technical measures:
 - Storage conditions (suitable, to be avoided): Avoid direct sunlight, high temperature, high humidity. Store in cool place (temperature: -20 ~ 35 degree C, humidity: 45 ~ 85%).
 - Incompatible products: Conductive materials, water, seawater, strong oxidizers and strong acids
 - Packing material (recommended, not suitable): Insulative and tear proof materials are recommended.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

- Engineering measures:
No engineering measure is necessary during normal use. In case of internal cell materials' leakage, operate the local exhaust or improve ventilation.
- Control parameters

Common chemical name / General name	ACGIH (2009)	
	TLV-TWA	BEI
Lithium transition metal oxidate	0.02mg/m ³ (as cobalt) * 0.2mg/m ³ (as manganese) * 0.2 mg/m ³ (as nickel) *	-
Aluminum	10mg/m ³ (metal coarse particulate) 5mg/m ³ (flammable powder) 5mg/m ³ (weld fume)	-
Carbon (Natural graphite) (Artificial graphite)	2mg/m ³ (inhalant coarse particulate)	-
Copper	0.2mg/m ³ (fume) 1.0mg/m ³ (a coarse particulate, Mist)	-
Organic electrolyte	-	-

ACGIH: American Conference of Governmental Industrial Hygienists, Inc.

TLV-TWA: Threshold Limit Value-Time Weighted Average concentration

BEI: Biological Exposure Indices

* Not every product includes all of these metals.

- Personal protective equipment
Respiratory protection: Respirator with air cylinder, dust mask
Hand protection: Protective gloves
Eye protection: Goggle or protective glasses designed to protect against liquid splashes
Skin and body protection: Working clothes with long sleeve and long trousers

9. PHYSICAL AND CHEMICAL PROPERTIES

- Appearance
Physical state: Solid
Form: Cylindrical or Prismatic or Prismatic (laminated)
Color: Metallic color or black(without tube if it has tube)
Odor: No odor
- pH: NA
- Specific temperatures/temperature ranges at which changes in physical state occur:
There is no useful information for the product as a mixture.
- Flash point: NA
- Explosion properties: NA
- Density: NA
- Solubility ,with indication of the solvent(s): Insoluble in water

10. STABILITY AND REACTIVITY

- Stability: Stable under normal use
- Hazardous reactions occurring under specific conditions
 - Conditions to avoid: When a battery cell is exposed to an external short-circuit, crushes, deformation, high temperature above 100 degree C, it will be the cause of heat generation and ignition. Direct sunlight and high humidity.
 - Materials to avoid: Conductive materials, water, seawater, strong oxidizers and strong acids.
 - Hazardous decomposition products: Acrid or harmful gas is emitted during fire.

- Chronic toxicity/Long term toxicity: By the long-term inhalation of coarse particulate or fume, it is possible to cause lung damage (aluminum lungs).

Graphite

- Acute toxicity: Unknown.
- Local effects: When it goes into one's eyes, it stimulates one's eyes; conjunctivitis, thickening of corneal epithelium or edematous inflammation palpebra may be caused.
- Chronic toxicity/Long term toxicity:
Since the long-term inhalation of high levels of graphite coarse particulate may become a cause of a lung disease or a tracheal disease.
- Carcinogenicity:
Graphite is not recognized as a cause of cancer by research organizations and natural toxic substance research organizations of cancer.

Copper

- Acute toxicity:
60-100mg sized coarse particulate causes a gastrointestinal disturbance with nausea and inflammation.
TDLo, hypodermic - Rabbit 375mg/kg
- Local effects:
Coarse particulate stimulates a nose and a tracheal.
When it goes into one's eyes, the symptom of the reddening and the pain is caused.
- Sensitization: Sensitization of the skin may be caused by long-term or repetitive contact.

Organic Electrolyte

- Acute toxicity:
LD₅₀, oral - Rat 2,000mg/kg or more
 - Local effects: Unknown.
 - Skin irritation study: Rabbit - Mild
 - eye irritation study: Rabbit - Very severe
-

12. ECOLOGICAL INFORMATION

- Persistence/degradability:

Since a battery cell and the internal materials remain in the environment, do not bury or throw out into the environment.

13. DISPOSAL CONSIDERATIONS

- Recommended methods for safe and environmentally preferred disposal:

Product (waste from residues)

Specified collection or disposal of lithium ion battery is required by the law like as "battery control law" in several nations. Collection or recycle of the battery is mainly imposed on battery's manufacturer or importer in the nations recycle is required.

Contaminated packaging

Neither a container nor packing is contaminated during normal use. When internal materials leaked from a battery cell contaminates, dispose as industrial wastes subject to special control.

14. TRANSPORT INFORMATION

In the case of transportation, avoid exposure to high temperature and prevent the formation of any condensation. Take in a cargo of them without falling, dropping and breakage. Prevent collapse of cargo piles and wet by rain. The container must be handled carefully. Do not give shocks that result in a mark of hitting on a cell. Please refer to Section 7-HANDLING AND STORAGE also.

UN regulation

- ID number: 3480
- Proper shipping name:
Lithium ion batteries
- Class: 9 *
- Packing group: II *

** However this product is defined as above, it is **not** recognized as "DANGEROUS GOODS" when its transport condition accords with instructions or provisions depend on region and transportation mode. About the instructions or provisions, please see descriptions in box brackets of following regulations.*

Regulation depends on region and transportation mode

- Worldwide, air transportation:
IATA-DGR ["packing instruction 965 section II"]
(When batteries are packaged with equipments or contained in equipments, refer packing instruction 966 or 967 instead of 965.)
 - Worldwide, sea transportation:
IMO-IMDG Code [special provision 188]
 - Europe, road transportation:
ADR [special provision 188]
-

15. REGULATORY INFORMATION

- Regulations specifically applicable to the product:

Wastes Disposal and Public Cleaning Law [Japan]
Law for Promotion of Effective Utilization of resources [Japan]
US Department of Transportation 49 Code of Federal Regulations [USA]

** About overlapping regulations, please refer to Section 14-TRANSPORT INFORMATION.*

16. OTHER INFORMATION

- This safety data sheet is offered an agency who handles this product to handle it safely.
- The agency should utilize this safety data sheet effectively (put it up, educate person in charge) and take proper measures.
- ***The information contained in this Safety data sheet is based on the present state of knowledge and current legislation.***
- This safety data sheet provides guidance on health, safety and environmental aspects of the product and should not be construed as any guarantee of technical performance or suitability for particular applications.

Reference

Chemical substances information: Japan Advanced Information center of Safety and Health
International Chemical Safety Cards (ICSCs): International Occupational Safety and Health Information Centre (CIS)
Dangerous Goods Regulations – 53rd Edition Effective 1 January 2012: International Air Transport Association (IATA)
IMDG Code - 2010 Edition: International Maritime Organization (IMO)
The European Agreement concerning the International Carriage of Dangerous Goods by Road – 2011:
The United Nations Economic Commission for Europe (UNECE)
RTECS (CD-ROM)
MSDS of raw materials prepared by the manufactures

First edition: Apr. 28, 2010

Prepared and approved by
Battery System Management Department
Lithium-Ion Battery Business Unit
Energy Company
Sanyo Electric Co., Ltd.

Test Verification Of Conformity

On the basis of the tests undertaken, the sample(s) of the below product have been found to comply with the essential requirements of the referenced specifications at the time the tests were carried out.

Verification/Report Number(s):	S1210038
Applicant Name/Address:	HEARTWAY MEDICAL PRODUCTS Co., Ltd. No.6 Gongyequ 25th Rd., Nantun Dist., Taichung City 408, Taiwan (R.O.C.)
Manufacture Name/Address:	RUBICON TECH Co., Ltd. 2F-12, No.210, Gongyequ 38th Rd., Xitun Dist, Taichung, Taiwan (R.O.C.)
Product(s) Tested:	Li-Ion Rechargeable Battery Pack
Ratings and principal characteristics:	24Vdc, 11.5Ah
Compliance with standards:	UN38.3
Model(s):	RT-705-1W
Brand name:	N/A
Issue Date:	2012. 12. XX

This Verification is for the exclusive use of AnTek's Client and is provided pursuant to the agreement between AnTek and its Client. AnTek's responsibility and liability are limited to the terms and conditions of the agreement. AnTek assumes no liability to any party, other than to the Client in accordance with the agreement, for any loss, expense or damage occasioned by the use of this Verification. Only the Client is authorized to copy or distribute this Verification. Any use of the AnTek name or one of its marks for the sale or advertisement of the tested material, product or service must first be approved in writing by AnTek. The observations and test results referenced from this Verification are relevant only to the sample tested. This Verification by itself does not imply that the material, product, or service is or has ever been under a AnTek certification program. The conclusion of the consignment test is only valid for the provided sample. In the test report, IATA DGR is referred to "Dangerous Goods Regulations" of International Air Transport Association; IMO IMDG Code is referred to "INTER- NATIONAL MARITIME DANGEROUS GOODS CODE" of International Maritime Organization (IMO); GB12268 is referred to National Standard named "List of dangerous goods" .

Joe Y. L. Lin
Manager
AnTek Certification Inc.
2012.12.XX

LITHIUM ION BATTERY SAFETY TESTING REPORT

Applicant:	HEARTWAY MEDICAL PRODUCTS Co., Ltd. No.6 Gongyequ 25th Rd., Nantun Dist., Taichung City 408, Taiwan (R.O.C.)
Product:	Li-Ion Rechargeable Battery Pack
Model:	RT-705-1W
Rating:	24Vdc, 11.5Ah
Test method & Criterion	United Nations "Recommendations on the TRANSPORT OF DANGEROUS GOODS"
Appearance	--
Verification Issuing Office Name	AnTek Certification Inc. 7F., No. 351, Yangguang St., Neihu District, Taipei City, Taiwan
Test Performed date:	Nov. 14, 2012– Dec. 03, 2012
Test Items:	Altitude Simulation, Thermal Test, Vibration, Shock, External Short Circuit, Overcharge
Conclusion:	The sample has passed the test items of UN38.3

Prepared by:

Reviewed by:

Gary Cheng


Rick Lin

Gary Cheng
Engineer

Rick Lin
Project engineer

Test Items

No	Name of Test Items	Standard requirement or the Clause Number of Standard	Test Result	Conclusion	Remark
1.	Altitude Simulation	UN Manual of Tests and Criteria Section (4th) 38.3 Test T.1	See Appendix 1	Passed	--
2.	Thermal Test	UN Manual of Tests and Criteria Section (4th) 38.3 Test T.2	See Appendix 2	Passed	--
3.	Vibration	UN Manual of Tests and Criteria Section (4th) 38.3 Test T.3	See Appendix 3	Passed	--
4.	Shock	UN Manual of Tests and Criteria Section (4th) 38.3 Test T.4	See Appendix 4	Passed	--
5.	External Short Circuit	UN Manual of Tests and Criteria Section (4th) 38.3 Test T.5	See Appendix 5	Passed	--
6.	Impact	UN Manual of Tests and Criteria Section (4th) 38.3 Test T.6	See Appendix 6	N/A	Only test for cell level
7.	Overcharge	UN Manual of Tests and Criteria Section (4th) 38.3 Test T.7	See Appendix 7	Passed	--
8.	Forced Discharge	UN Manual of Tests and Criteria Section (4th) 38.3 Test T.8	See Appendix 8	N/A	Only test for cell level
Test Environment Condition		Ambient Temperature: 23.8 °C Ambient Humidity: 54.6%			

*Test 6 and Test 8 are required to cell level, not battery pack level.

Test Procedure

Pre-Learning:

- (a) When testing primary cells and batteries under tests 1 to 5 the following shall be tested in the quantity indicated:
- (i) ten cells in undischarged states;
 - (ii) ten cells in fully discharged states;
 - (iii) four small batteries in undischarged states;
 - (iv) four small batteries in fully discharged states;
 - (v) four large batteries in undischarged states; and
 - (vi) four large batteries in fully discharged states.
- (b) When testing rechargeable cells and batteries under tests 1 to 5 the following shall be tested in the quantity indicated:
- (i) ten cells at first cycle, in fully charged states;
 - (ii) four small batteries at first cycle, in fully charged states;
 - (iii) four small batteries after 50 cycles ending in fully charged states;
 - (iv) two large batteries at first cycle, in fully charged states; and
 - (v) two large batteries after 25 cycles ending in fully charged states.
- (c) When testing primary and rechargeable cells under test 6, the following shall be tested in the quantity indicated:
- (i) for primary cells, five cells in undischarged states and five cells in fully discharged states;
 - (ii) for component cells of primary batteries, five cells in undischarged states and five cells in fully discharged states;
 - (iii) for rechargeable cells, five cells at first cycle at 50% of the design rated capacity; and
 - (iv) for component cells of rechargeable batteries, five cells at first cycle at 50% of the design rated capacity.
- For prismatic cells, ten test cells are required instead of the five described above, so that the procedure can be carried out on five cells along the longitudinal axes and, separately, five cells along the other axes. In every case, the test cell is only subjected to one impact.
- (d) When testing rechargeable batteries under test 7, the following shall be tested in the quantity indicated:
- (i) four small batteries at first cycle, in fully charged states;
 - (ii) four small batteries after 50 cycles ending in fully charged states;
 - (iii) two large batteries at first cycle, in fully charged states; and
 - (iv) two large batteries after 25 cycles ending in fully charged states.
- Batteries not equipped with overcharge protection that are designed for use only in a battery assembly, which affords such protection, are not subject to the requirements of this test.

- (e) When testing primary and rechargeable cells under test 8, the following shall be tested in the quantity indicated:
- (i) ten primary cells in fully discharged states;
 - (ii) ten rechargeable cells, at first cycle in fully discharged states; and
 - (iii) ten rechargeable cells after 50 cycles ending in fully discharged states.
- (f) When testing a battery assembly in which the aggregate lithium content of all anodes, when fully charged, is not more than 500 g, or in the case of a lithium ion battery, with a Watt-hour rating of not more than 6 200 Watt-hours, that is assembled from cells or batteries that have passed all applicable tests, one battery assembly in a fully charged state shall be tested under tests 3, 4 and 5, and, in addition, test 7 in the case of a rechargeable battery assembly. For a rechargeable battery assembly, the assembly shall have been cycled at least 25 cycles.

Test 1: Altitude simulation test

Test procedure:

Test cells and batteries shall be stored at a pressure of 11.6 kPa or less for at least six hours at ambient temperature (20 ± 5 °C).

Requirement:

Cells and batteries meet this requirement if there is no mass loss, no leakage, no venting, no disassembly, no rupture and no fire and if the open circuit voltage of each test cell or battery after testing is not less than 90% of its voltage immediately prior to this procedure. The requirement relating to voltage is not applicable to test cells and batteries at fully discharged states.

Test 2: Thermal test

Test procedure:

Test cells and batteries are to be stored for at least six hours at a test temperature equal to 75 ± 2 °C, followed by storage for at least six hours at a test temperature equal to -40 ± 2 °C. The maximum time interval between test temperature extremes is 30 minutes. This procedure is to be repeated 10 times, after which all test cells and batteries are to be stored for 24 hours at ambient temperature (20 ± 5 °C). For large cells and batteries the duration of exposure to the test temperature extremes should be at least 12 hours.

Requirement:

Cells and batteries meet this requirement if there is no mass loss, no leakage, no venting, no disassembly, no rupture and no fire and if the open circuit voltage of each test cell or battery after testing is not less than 90% of its voltage immediately prior to this procedure. The requirement relating to voltage is not applicable to test cells and batteries at fully discharged states.

Test 3: Vibration test**Test procedure:**

Cells and batteries are firmly secured to the platform of the vibration machine without distorting the cells in such a manner as to faithfully transmit the vibration. The vibration shall be a sinusoidal waveform with a logarithmic sweep between 7 Hz and 200 Hz and back to 7 Hz traversed in 15 minutes. This cycle shall be repeated 12 times for a total of 3 hours for each of three mutually perpendicular mounting positions of the cell. One of the directions of vibration must be perpendicular to the terminal face.

The logarithmic frequency sweep is as follows: from 7 Hz a peak acceleration of 1 gn is maintained until 18 Hz is reached. The amplitude is then maintained at 0.8 mm (1.6 mm total excursion) and the frequency increased until a peak acceleration of 8 gn occurs (approximately 50 Hz). A peak acceleration of 8 gn is then maintained until the frequency is increased to 200 Hz.

Requirement:

Cells and batteries meet this requirement if there is no mass loss, no leakage, no venting, no disassembly, no rupture and no fire and if the open circuit voltage of each test cell or battery after testing is not less than 90% of its voltage immediately prior to this procedure. The requirement relating to voltage is not applicable to test cells and batteries at fully discharged states.

Test 4: Shock test**Test procedure:**

Test cells and batteries shall be secured to the testing machine by means of a rigid mount which will support all mounting surfaces of each test battery. Each cell or battery shall be subjected to a half-sine shock of peak acceleration of 150 gn and pulse duration of 6 milliseconds. Each cell or battery shall be subjected to three shocks in the positive direction followed by three shocks in the negative direction of three mutually perpendicular mounting positions of the cell or battery for a total of 18 shocks.

However, large cells and large batteries shall be subjected to a half-sine shock of peak acceleration of 50 gn and pulse duration of 11 milliseconds. Each cell or battery is subjected to three shocks in the positive direction followed by three shocks in the negative direction of each of three mutually perpendicular mounting positions of the cell for a total of 18 shocks.

Requirement:

Cells and batteries meet this requirement if there is no mass loss, no leakage, no venting, no disassembly, no rupture and no fire and if the open circuit voltage of each test cell or battery after testing is not less than 90% of its voltage immediately prior to this procedure. The requirement relating to voltage is not applicable to test cells and batteries at fully discharged states.

Test 5: External short circuit test**Test procedure:**

The cell or battery to be tested shall be temperature stabilized so that its external case temperature reaches 55 ± 2 °C and then the cell or battery shall be subjected to a short circuit condition with a total external resistance of less than 0.1 ohm at 55 ± 2 °C. This short circuit condition is continued for at least one hour after the cell or battery external case temperature has returned to 55 ± 2 °C. The cell or battery must be observed for a further six hours for the test to be concluded.

Requirement:

Cells and batteries meet this requirement if their external temperature does not exceed 170 °C and there is no disassembly, no rupture and no fire within six hours of this test.

Test 6: Impact test**Test procedure:**

The test sample cell or component cell is to be placed on a flat surface. A 15.8 mm diameter bar is to be placed across the centre of the sample. A 9.1 kg mass is to be dropped from a height of 61 ± 2.5 cm onto the sample.

A cylindrical or prismatic cell is to be impacted with its longitudinal axis parallel to the flat surface and perpendicular to the longitudinal axis of the 15.8 mm diameter curved surface lying across the centre of the test sample. A prismatic cell is also to be rotated 90 degrees around its longitudinal axis so that both the wide and narrow sides will be subjected to the impact. Each sample is to be subjected to only a single impact. Separate samples are to be used for each impact.

A coin or button cell is to be impacted with the flat surface of the sample parallel to the flat surface and the 15.8 mm diameter curved surface lying across its centre.

Requirement:

Cells and component cells meet this requirement if their external temperature does not exceed 170 °C and there is no disassembly and no fire within six hours of this test.

Test 7: Overcharge test**Test procedure:**

The charge current shall be twice the manufacturer's recommended maximum continuous charge current. The minimum voltage of the test shall be as follows:

- (a) when the manufacturer's recommended charge voltage is not more than 18V, the minimum voltage of the test shall be the lesser of two times the maximum charge voltage of the battery or 22V.
- (b) when the manufacturer's recommended charge voltage is more than 18V, the minimum