# SPORTleg

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# Manual

**Español**  
04

ES

**Portugues**  
12

PO

**English**  
20

EN

**Français**  
28

FR

**Italiano**  
36

IT

**Deutsch**  
44

DE

Estimado consumidor,

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Dear customer,


Thank you for acquiring the **Vida10 SPORTleg** training system.

For already more than 15 years, helping you to take care of your wellbeing is our main and most important aim. We provide you therefore with yet another personal domestic training device, conceived thanks to the latest technology, easy and practical to use, so that you may treat you to self-care wherever you are, whenever you need it.

Please store this manual in a safe place for future reference.

Please remember:

1. The photos of the product correspond to the latest model and are for reference purposes only.
2. This manual contains recommendations on how to use this product.

#### Correct disposal of this product

This symbol indicates that, within the EU, this product should not be disposed of along with other household waste. To prevent potential damage to the environment or human health from the uncontrolled elimination of waste, please recycle it responsibly to promote the sustainable reuse of its material resources. To return your used product, use the return and collection channels or contact the seller. They can take the product to the recycling point.

# index

- 22** Introduction
- 23** How to use the leg exerciser
- 24** Remote control functions
- 24** Equivalence of each level
- 25** Usage recommendations
- 25** Exercise on a daily
- 27** Datos técnicos

## **Do you spend hour after hour, day after day sitting or standing without doing any kind of exercise?**

This lack of movement may afflict your daily wellbeing, disavouring your positive attitude and feeling of fulfillment.

If you spend long hours each day sitting in the same chair at home or at work, you will have noticed that your legs are weakened and feel uncomfortable after having sat in the same position for a prolonged time.

The personal leg exerciser Sportleg is the only motorized solution of a passive exerciser that creates a constant leg movement to stimulate your lower extremities while you remain seated.

It allows you to comfortably train your legs, while traditional exercise techniques require space and positions that, due to various difficulties, for many people are impossible to realise.

## How to use the leg exerciser

1. Connect the power adapter with the device and plug it into a power socket.

2. Press the ON/OFF button (B) to turn on the device; press the + or – button to increase or decrease the speed between level 1 and level 5.

3. Sit down and put your feet on the pedals (A).

4. Having finished your training, turn off the device.

### **Optional**

To incline the device upwards, pull out the legs (C) from the underside.

## Remote control functions

## Equivalence of each level

**Level 1:** equals 32 steps per minute.

**Level 2:** equals 36 steps per minute.

**Level 3:** equals 40 steps per minute.

**Level 4:** equals 44 steps per minute.

**Level 5:** equals 52 steps per minute.

## Usage recommendations

- Do not sit or stand on the device.
- Keep the device out of reach of children.
- This device does not have any playful purpose and must not be used as a toy.
- Do not use the device on a slippery or inclined surface.
- Do not wet the device.
- Do not put fingers or toes between the pedals.

## Exercise on a daily basis to improve your wellbeing

**Sportleg** has been especially designed to improve the wellbeing of people who move very little in their daily lives and have to remain seated for a long time each day.

A constant use in moments when we have to stay seated while performing a parallel task may be the best occasion to use **Sportleg** and benefit from its positive effect on our sensation of wellbeing.

More and more professionals find themselves forced to remain seated during their entire working day. The leg exerciser is the best tool to prevent discomfort and pain caused by a lack of exercise.

Thanks to this device, your legs stay active because of the constant leg movement it creates. This movement visibly improves the positive feeling of our body.

In addition to that, since device is a passive exercising system, you will train your legs without any effort. For example, 3 hours of exercise have a similar result as 1 hour of walking.

Keeps your body active

Helps to stimulate your legs after sitting for a long time

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