International patents pending. All imitations will be prosecuted.
Gymform® Leg Action is a registered EU/CTM trade mark

-  **Workout Guide**
-  **Guía de ejercicios**
-  **Guide de séance d'entraînement**
-  **Training Führer**
-  **Guida di allenamento**
-  **Guia do exercicio**
-  **Training gids**

English

Note: If doing a particular massage is not comfortable, stop the massage. Use Gymform® Leg Action only if it feels comfortable. For massages, we recommend always using programs that provide both massage and vibration. Choose a program with vibration and massage (P2, P4 or M1).

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English

Time of use: Maximum 10 minutes.

Benefits: This massage helps relax the feet and legs while activating circulation and relieving muscle tension.

Position: Sit comfortably in a chair or on the sofa with your back supported and place your feet on the foot pads

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
English

Time of use: Maximum 10 minutes.
Benefits: This massage is designed to work on the calves and gently massage the feet and thighs.

Position: Sit with your back straight on the edge of a chair or sofa. Place the device as far away as possible while still keeping the soles of your feet completely on the foot pads. Your calves should be stretched.

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English

Time of use: Maximum 10 minutes.

Benefits: This massage is designed to work the soles of the feet, calves and joints. It is a strong, deep massage.

Position: Sit with your back straight on the edge of a chair or the sofa. Move the device closer to your body and place your feet on the foot pads.

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English

Time of use: 10 minutes

Benefits: This massage strengthens the calves.

Position: Sit in a chair or on the sofa with your back straight.

Step 1: Place your toes on the bottom part of the foot pad and keep your heels at the same level as your toes.

Step 2: Raise your heels and tense your calves. Hold this position for 5 seconds.

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English

Time of use: 10 minutes

Benefits: This massage helps relax the calves.

Position: Sit on the floor and place your calves on the foot pads.

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English

Time of use: Maximum 10 minutes.

Benefits: This massage is designed to work the hamstring muscles.

Position: Sit on the floor and place your thighs on the foot pads. We recommend using Manual Programme M1 to choose the level of intensity of the massage you've chosen.

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English

Time of use: Maximum 10 minutes.

Benefits: This massage helps relieve metatarsal pain.

Position: Sit in a comfortable position and place the painful area on top of the shiatsu massage rollers. Exert moderate pressure on the area. We recommend using Manual Programme M1 at slow speeds. seconds.

Rest for 5 seconds and then do 10 (beginners), 15 (advanced users) or 20 repetitions (experienced users).

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
English

Time of use: Maximum 10 minutes.
Benefits: This massage helps relax the thighs and relieve thigh pain.
Position: Assume the position shown in the picture and then do the same massage on the other leg.
Use Manual Programme M1 to adjust to the exact speed of your choice.
Important: Do not let your entire body weight fall on the device. It may break.

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English

Time of use: Maximum 10 minutes.
Benefits: This massage helps relax the thighs and relieve thigh pain.

Position: Assume the position shown in the picture. Use Manual Programme M1 to adjust to the exact speed of your choice.

Important: Do not let your entire body weight fall on the device. It may break.

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Lower Back Massage ***Masaje en zona lumbar*** ***Massage lombaire***

English

Time of use: Maximum 10 minutes.
Benefits: This massage helps relax the lower back and relieve lower back pain.

Position: Position the Gymform® Leg Action against the chair back so the shiatsu massage rollers are in contact with the lower back. Lean back against the device so you are comfortable and ready for the massage.

Use Manual Programme M1 to adjust to the exact speed of your choice.

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ISL GmbH
Best Direct (International) Benelux B.V.
Best Direct International Ltd,
ISL ITALY S.r.l.